

Anti-TAOK1 Rabbit antibody



Catalog Number: 167627

Product name

Anti-TAOK1 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IP

Antibody description

TAOK1 Rabbit monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic peptide of human TAOK1

IP: 1/200

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of TAOK1 in C6,Ramos
using TAOK1 Rabbit mAb(1:1000 diluted)