

CCT3 antibody



Catalog Number: 107037

Product name

CCT3 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

CCT3 Mouse Monoclonal antibody. Positive WB detected in COLO 320 cells. Positive IF detected in HepG2 cells. Positive IHC detected in human testis tissue, human liver cancer tissue. Positive FC detected in HepG2 cells. Observed molecular weight by Western-blot: 56-60kd

Preparation

This antibody was obtained by immunization of CCT3 recombinant protein (Accession Number: NM_005998). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB, IF, IHC, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:107037(CCT3 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:107037 (CCT3 Antibody) at dilution of 1:50 and Rhodamine-labeled goat anti-mouse IgG




(red).

Immunohistochemistry of paraffin-embedded


CCT3 antibody

Catalog Number: 107037

human testis tissue slide using Catalog
No:107037(CCT3 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded
human testis tissue slide using Catalog
No:107037(CCT3 Antibody) at dilution of 1:50



1X10⁶ HepG2 cells were stained with 0.2ug
CCT3 antibody (Catalog No:107037, red) and
control antibody (blue). Fixed with 90% MeOH
blocked with 3% BSA (30 min). Alexa Fluor 488-
conjugated AffiniPure Goat Anti-Mouse IgG(H+L)
with dilution 1:1000.