

TSNAX antibody



Catalog Number: 116437

Product name

TSNAX antibody

Immunogen

[Human TSNAX Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TSNAX Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IF detected in MCF-7 cells.

Preparation

This antibody was obtained by immunization of TSNAX recombinant protein (Accession Number: NM_005999). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116437(TSNAX antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of MCF-7 cells, using TSNAX antibody Catalog No:116437 at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG (green). Blue pseudocolor = DAPI (fluorescent DNA dye).