

ORM1/2 antibody



Catalog Number: 107437

Product name

ORM1/2 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human; other species not tested.

Antibody description

ORM1/2 Mouse Monoclonal antibody. Positive IHC detected in human liver tissue, human liver cancer tissue. Positive IF detected in HepG2 cells. Positive WB detected in Human Blood. Observed molecular weight by Western-blot: 41kDa

Preparation

This antibody was obtained by immunization of ORM1/2 recombinant protein (Accession Number: NM_000607). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:107437(ORM1 antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:107437(ORM1 Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells, using ORM1 antibody Catalog No: at 1:25 dilution and

ORM1/2 antibody



Catalog Number: 107437

Rhodamine-labeled goat anti-mouse IgG (red).