

DYNC1LI1 antibody



Catalog Number: 110137

Product name

DYNC1LI1 antibody

Immunogen

[Human DYNC1LI1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

DYNC1LI1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in K-562 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 57-60 kDa

Preparation

This antibody was obtained by immunization of DYNC1LI1 recombinant protein (Accession Number: NM_016141). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110137(DYNC1LI1 Antibody) at dilution of



1:1000

Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:110137(DYNC1LI1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)