

Anti-p44/p42 MAPK (ERK1/2) (4F3) Mouse antibody



Catalog Number: 168137

Product name

Monoclonal

Anti-p44/p42 MAPK (ERK1/2) (4F3) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

p44/p42 MAPK (ERK1/2) (4F3) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Human p44 MAPK (Erk1) synthetic peptide conjugated to KLH.

Formulation


Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Validations



Western blot analysis of extracts from Jurkat, A549, C6, 3T3 and HeLa cell lysates using p44/42 MAPK (Erk1/2) mouse mAb (1:1000 diluted). Predicted band size: 42/44 kDa. Observed band size: 42/44 kDa.