

RANBP10 antibody



Catalog Number: 114537

Product name

RANBP10 antibody

IHC: 1:20-1:200

Specificity

human, mouse, rat; other species not tested.

Antibody description

RANBP10 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human skeletal muscle tissue, MCF7 cells, mouse heart tissue. Positive IHC detected in human hepatocirrhosis. Observed molecular weight by Western-blot: 67-70kd

Preparation

This antibody was obtained by immunization of RANBP10 recombinant protein (Accession Number: NM_020850). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114537(RANBP10 antibody) at dilution of



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:114537(RANBP10 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog

RANBP10 antibody



Catalog Number: 114537

No:114537(RANBP10 antibody) at dilution of 1:50 (under 40x lens)