

Anti-KDM1/LSD1 Rabbit antibody



Catalog Number: 166937

Product name

Anti-KDM1/LSD1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

KDM1/LSD1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human KDM1/LSD1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/20-1/100


FC: 1/20

IP: 1/20-1/50

Validations




Immunocytochemistry of KDM1/LSD1(green) in Hela cells using KDM1/LSD1 Rabbit pAb at



dilution 1/50, and DAPI(blue)

Western blot detection of KDM1/LSD1 in Hela, CHO-K1, C6, Jurkat cell lysates using KDM1/LSD1 Rabbit pAb(1:500 diluted). Predicted band size: 93kDa. Observed band size: 110kDa.



Western blot detection of KDM1/LSD1 in Jurkat cell lysates using KDM1/LSD1 Rabbit pAb(1:1000 diluted). Predicted band size: 93KDa. Observed band size: 110KDa.