

# NKX3-2 antibody



Catalog Number: 113237

---

## Product name

NKX3-2 antibody

## Immunogen

[Human NKX3-2 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

NKX3-2 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 36-40 kDa

## Preparation

This antibody was obtained by immunization of NKX3-2 recombinant protein (Accession Number: NM\_001189). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113237(NKX3-2 Antibody) at dilution of 1:300