

Product name

TRAIP antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

IF: 1:10-1:100

Antibody description

TRAIP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A549 cells, HepG2 cells, human brain tissue, human placenta tissue, mouse placenta tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of TRAIP recombinant protein (Accession Number: NM_005879). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions


Recommended Dilution:

Validations

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116237(TRAIP antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:116237(TRAIP Antibody) at dilution of 1:200 (under 10x lens).




TRAIP antibody



Catalog Number: 116237

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:116237(TRAIP Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of HepG2 cells using Catalog No:116237(TRAIP Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)