

Anti-VCP Rabbit antibody



Catalog Number: 167237

Product name

Anti-VCP Rabbit antibody

IHC: 1/20-1/100

ICC/IF: 1/20-1/50

Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20-1/50

Antibody description

VCP Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human VCP

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/2000-1/10000



Western blot detection of VCP in MCF7, C6, 3T3 cell lysates using VCP Rabbit pAb(1:1000 diluted). Predicted band size: 89KDa. Observed band size: 89KDa.



Immunohistochemistry of VCP in paraffin-embedded Human tonsil using VCP Rabbit pAb at dilution 1/20