

# TET3 Polyclonal Antibody



Catalog Number: 161637

---

## Product name

TET3 Polyclonal Antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

IHC, IF

## Antibody description

Polyclonal antibody to TET3

## Dilutions


IHC 1:50 - 1:200

IF 1:50 - 1:200

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1556-1795 of human TET3 (NP\_001274420.1).

## Validations



## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

Immunofluorescence - TET3 Polyclonal Antibody

## Clonality

Polyclonal

Immunofluorescence analysis of A549 cells using TET3 antibody . Blue: DAPI for nuclear staining.