

BCKDK Polyclonal Antibody



Catalog Number: 165737

Product name

BCKDK Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BCKDK

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 153-412 of human BCKDK (NP_005872.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations




Western blot - BCKDK Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using BCKDK antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - BCKDK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using BCKDK antibody at dilution of 1:100



(40x lens).


BCKDK Polyclonal Antibody



Catalog Number: 165737

Immunohistochemistry - BCKDK Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using BCKDK antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - BCKDK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using BCKDK antibody at



dilution of 1:100 (40x lens).

Immunofluorescence - BCKDK Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using BCKDK antibody at dilution of 1:100. Blue: DAPI for nuclear staining.