

Anti-SOX9 Rabbit antibody



Catalog Number: 168737

Product name

Anti-SOX9 Rabbit antibody

Applications

WB, IHC-P, ICC/IF, FC, IP

Specificity

Human, Mouse, Rat

Dilutions

1:2000

Antibody description

1:200

SOX9 Rabbit monoclonal antibody

1:50

Preparation

1:20

Antigen: Recombinant protein of human SOX9

1:20

Formulation

Validations

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG



Western blot detection of SOX9 in 3T3,HeLa cell lysates using SOX9 Rabbit mAb(1:1000 diluted).Predicted band size:56kDa.Observed band size:70kDa.