

ZFYVE27 antibody



Catalog Number: 117047

Product name

ZFYVE27 antibody

Immunogen

[Human ZFYVE27 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ZFYVE27 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse thymus tissue, Jurkat cells, mouse brain tissue, mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of ZFYVE27 recombinant protein (Accession Number: NM_144588). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:117047(ZFYVE27 antibody) at dilution of



1:500

IP Result of anti-ZFYVE27 (IP:Catalog No:117047, 3ug; Detection:Catalog No:117047 1:500) with mouse testis tissue lysate 4000ug.