

General Information

Protein Construction

A DNA sequence encoding the human BIRC5 (NP_001159.2) (Met1-Asp142) was expressed and purified with two additional amino acids (Gly & Pro) at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly

Molecular Mass

The recombinant human BIRC5 consists of 144 amino acids and predicts a molecular mass of 16.5 KDa. It migrates as an approximately 19 KDa band

in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, 20% glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human Survivin / BIRC5 / API4 Protein SDS-PAGE