

SPATA1 antibody



Catalog Number: 115547

Product name

SPATA1 antibody

Immunogen

[Human SPATA1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

SPATA1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 40-50 kDa

Preparation

This antibody was obtained by immunization of SPATA1 recombinant protein (Accession Number: BC127634). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115547(SPATA1 Antibody) at dilution of 1:300