

TERF2 antibody



Catalog Number: 115947

Product name

TERF2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TERF2 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in Jurkat cells, HeLa cells, K-562 cells, MCF7 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 65-69 kDa

Preparation

This antibody was obtained by immunization of TERF2 recombinant protein (Accession Number: BC024890). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115947(TERF2 antibody) at dilution of 1:1000



IP Result of anti-TERF2 (IP:Catalog No:115947, 4ug; Detection:Catalog No:115947 1:500) with K-562 cells lysate 4000ug.