

# MPP1 antibody



Catalog Number: 112747

---

**Product name**

pH 7.3.

MPP1 antibody

**Storage****Immunogen**

Store at -20°C. DO NOT ALIQUOT

[Human MPP1 Recombinant protein \(GST tag & His tag\)](#)

**Clonality**

Polyclonal

**Specificity**

human, mouse, rat; other species not tested.

**Ig Type**

Rabbit IgG

**Antibody description**

MPP1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, mouse thymus tissue. Observed molecular weight by Western-blot: 55kd

**Applications**

ELISA, WB

**Preparation**

This antibody was obtained by immunization of MPP1 recombinant protein (Accession Number: NM\_002436). Purification method: Antigen affinity purified.

**Dilutions**

Recommended Dilution:  
WB: 1:500-1:5000

**Formulation**

PBS with 0.02% sodium azide and 50% glycerol

**Validations**

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112747(MPP1 antibody) at dilution of 1:600