

# Anti-CDKN2A/p16INK4a Rabbit antibody



Catalog Number: 166747

## Product name

Rabbit IgG

Anti-CDKN2A/p16INK4a Rabbit antibody

## Applications

## Specificity

WB, FC, IP

Human

## Dilutions

## Antibody description

WB: 1/1000

CDKN2A/p16INK4a Rabbit polyclonal antibody

FC: 1/20

## Preparation

IP: 1/20

Antigen: Recombinant protein of human CDKN2A/p16INK4a

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type



Western blot detection of CDKN2A/p16INK4a in 3T3,HeLa cell lysates using CDKN2A/p16INK4a Rabbit pAb(1:1000 diluted).Predicted band size:17kDa.Observed band size:17kDa.