

RAB6B antibody



Catalog Number: 114457

Product name

RAB6B antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

RAB6B Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, C6 cells, rat brain tissue. Positive IF detected in C6 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of RAB6B recombinant protein (Accession Number: NM_016577). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114457(RAB6B antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114457(RAB6B antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded

RAB6B antibody


Catalog Number: 114457

human gliomas using Catalog No:114457(RAB6B antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-RAB6B (IP:Catalog No:114457, 3ug; Detection:Catalog No:114457 1:1000) with

mouse brain tissue lysate 4000ug.



Immunofluorescent analysis of C6 cells using Catalog No:114457(RAB6B Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)