

MRPS34 antibody



Catalog Number: 112857

Product name

MRPS34 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MRPS34 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human heart tissue. Positive IHC detected in human heart tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of MRPS34 recombinant protein (Accession Number: NM_023936). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112857(MRPS34 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:112857(MRPS34 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human heart using Catalog No:112857(MRPS34 antibody) at dilution of 1:50 (under 10x lens)

MRPS34 antibody



Catalog Number: 112857

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of Hela cells, using MRPS34 antibody Catalog No:112857 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).