

Grancalcin antibody



Catalog Number: 111157

Product name

Grancalcin antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Grancalcin Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Positive IHC detected in human tonsillitis tissue, human spleen tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of Grancalcin recombinant protein (Accession Number: NM_012198). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:111157(GCA Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:111157(GCA Antibody) at dilution of 1:50



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:111157(GCA Antibody) at dilution of 1:600