

SFTPb antibody



Catalog Number: 115157

Product name

SFTPb antibody

Specificity

Human; other species not tested.

Antibody description

SFTPb Rabbit Polyclonal antibody. Positive IP detected in Y79 cells. Positive WB detected in Y79 cells. Positive IHC detected in human heart tissue. Observed molecular weight by Western-blot: 43kd

Preparation

This antibody was obtained by immunization of SFTPb recombinant protein (Accession Number: NM_198843). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:1000

IHC: 1:20-1:200


Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:115157(SFTPb antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:115157(SFTPb antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-SFTPb (IP:Catalog No:115157, 3ug; Detection:Catalog No:115157 1:300) with Y79 cells lysate 3000ug.