

S100A13 antibody



Catalog Number: 114957

Product name

S100A13 antibody

Recommended Dilution:

WB: 1:200-1:1000

Specificity


IHC: 1:20-1:200

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

S100A13 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells. Positive IHC detected in human thyroid cancer tissue, human brain tissue, human heart tissue, human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue, human thyroid tissue. Observed molecular weight by Western-blot: 11 kDa




Preparation

This antibody was obtained by immunization of S100A13 recombinant protein (Accession Number: NM_001024210). Purification method: Antigen affinity purified.

MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114957(S100A13 Antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemistry of paraffin-embedded human thyroid cancer tissue slide using Catalog No:114957(S100A13 Antibody) at dilution of 1:50 (under 10x lens)

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions



S100A13 antibody



Catalog Number: 114957

Immunohistochemistry of paraffin-embedded
human thyroid cancer tissue slide using Catalog

No:114957(S100A13 Antibody) at dilution of 1:50
(under 40x lens)