

Anti-SAE2/UBA2 Rabbit antibody



Catalog Number: 168957

Product name

Rabbit IgG

Anti-SAE2/UBA2 Rabbit antibody

Applications

Specificity

WB, ICC/IF, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

SAE2/UBA2 Rabbit monoclonal antibody

ICC/IF: 1/20

Preparation

IP: 1/20

Antigen: A synthetic peptide of human UBA2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Validations



Western blot detection of SAE2 / UBA2 in Rat Brain, C6, HeLa cell lysates using SAE2 / UBA2 Rabbit mAb(1:1000 diluted). Predicted band size: 71 kDa. Observed band size: 90 kDa.