

NMUR1 antibody



Catalog Number: 113257

Product name

NMUR1 antibody

IHC: 1:20-1:200


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

NMUR1 Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue, human colon tissue, human testis tissue. Positive WB detected in human colon tissue, human stomach tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 47kd




Preparation

This antibody was obtained by immunization of NMUR1 recombinant protein (Accession Number: NM_006056). Purification method: Antigen affinity purified.

human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:113257(NMUR1 antibody) at dilution of 1:300

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:113257(NMUR1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:113257(NMUR1 Antibody) at dilution of 1:50

NMUR1 antibody



Catalog Number: 113257

(under 40x lens)