

UBP1 antibody



Catalog Number: 116657

Product name

UBP1 antibody

Immunogen

[Human UBP1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

UBP1 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue. Observed molecular weight by Western-blot: 68kd

Preparation

This antibody was obtained by immunization of UBP1 recombinant protein (Accession Number: NM_014517). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116657(UBP1 antibody) at dilution of 1:1000