

NUDT16 antibody



Catalog Number: 113357

Product name

NUDT16 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

IF: 1:20-1:200

Antibody description

NUDT16 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive IHC detected in human kidney tissue. Positive WB detected in HeLa cells, HEK-293 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of NUDT16 recombinant protein (Accession Number: NM_152395). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP, IF

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113357(NUDT16 antibody) at dilution of 1:2500



IP Result of anti-NUDT16 (IP:Catalog No:113357, 3ug; Detection:Catalog No:113357 1:1000) with HEK-293 cells lysate 2400ug.




NUDT16 antibody




Catalog Number: 113357

Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:113357(NUDT16 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog

No:113357(NUDT16 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HEK-293 cells using Catalog No:113357(NUDT16 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)