

# Rolinsatamab reference antibody



Catalog Number: 177067

## Product name

Rolinsatamab reference antibody

## Specificity

Human PRLR / Prolactin Receptor

## Antibody description

Anti-PRLR / Prolactin Receptor Antibody (Rolinsatamab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications

ELISA

## Validations



## SDS-PAGE

Anti-PRLR / Prolactin Receptor Antibody (rolinsatamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.