

Anti-MUM1 Rabbit antibody



Catalog Number: 167867

Product name

Rabbit IgG

Anti-MUM1 Rabbit antibody

Applications

Specificity

WB, IHC-P, IP

Human

Dilutions

Antibody description

WB: 1/1000

MUM1 Rabbit monoclonal antibody

IHC: 1/50

Preparation

IP: 1/20

Antigen: A synthetic peptide of human MUM1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage


Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Validations



Western blot detection of MUM1 in Ramos cell lysates using MUM1 Rabbit mAb(1:1000 diluted).Predicted band size:52kDa.Observed band size:52kDa.