

# CDK8 antibody

Catalog Number: 109167

## Product name

CDK8 antibody

IHC: 1:20-1:200

IF: 1:50-1:500


## Specificity

Human, Mouse; other species not tested.

## Validations

## Antibody description

CDK8 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, HEK-293 cells, HeLa cells. Positive IF detected in NIH/3T3 cells. Positive IHC detected in human stomach cancer tissue, human testis tissue. Observed molecular weight by Western-blot: 53 kDa




## Preparation

This antibody was obtained by immunization of CDK8 recombinant protein (Accession Number: NM\_001318368). Purification method: Antigen affinity purified.

NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109167(CDK8 antibody) at dilution of 1:500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human stomach cancer using Catalog No:109167(CDK8 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IF



## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded


# CDK8 antibody



Catalog Number: 109167

---

human stomach cancer using Catalog No:109167(CDK8 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (10% Formaldehyde) fixed NIH/3T3 cells using Catalog No:109167(CDK8 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)