

Anti-Cystathionase/CTH Rabbit antibody



Catalog Number: 167567

Product name

Rabbit IgG

Anti-Cystathionase/CTH Rabbit antibody

Applications

Specificity

WB, IP

Human

Dilutions

Antibody description

WB: 1/2000-1/10000

Cystathionase/CTH Rabbit monoclonal antibody

IP: 1/20


Preparation

Antigen: A synthetic peptide of human Cystathionase/CTH

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Western blot detection of Cystathionase/CTH in HeLa,A549,HL-60 cell lysates using Cystathionase/CTH Rabbit mAb(1:1000 diluted).Predicted band size:45kDa.Observed band size:45kDa.