

ZSCAN21 antibody



Catalog Number: 117267

Product name

ZSCAN21 antibody

Immunogen

[Human ZSCAN21 Recombinant protein \(GST tag & His tag\)](#)

Specificity

human; other species not tested.

Antibody description

ZSCAN21 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 54kd

Preparation

This antibody was obtained by immunization of ZSCAN21 recombinant protein (Accession Number: NM_001362779). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117267(ZSCAN21 antibody) at dilution of 1:1000