

Anti-H-FABP (3E8B10) Mouse antibody



Catalog Number: 168267

Product name

Mouse IgG2b

Anti-H-FABP (3E8B10) Mouse antibody

Applications

Specificity

sELISA Detector

Human

Validations

Antibody description

H-FABP (3E8B10) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant H-FABP protein expressed in E.coli


Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Standard Curve for H-FABP: Capture Antibody Mouse mAb (5C4-F5) to H-FABP at 4u03bcg/ml and Detector Antibody Mouse mAb(3E8-B10u30017B10-G11u30014H1-B8)to H-FABP at 0.08 u03bcg/ml.