

QKI antibody



Catalog Number: 114367

Product name

QKI antibody

Immunogen

[Human QKI Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

QKI Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 38kd

Preparation

This antibody was obtained by immunization of QKI recombinant protein (Accession Number: NM_006775). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114367(QKI antibody) at dilution of 1:400



IP Result of anti-QKI (IP:Catalog No:114367, 4ug; Detection:Catalog No:114367 1:400) with K-562 cells lysate 3600ug.