

# TSPY3 antibody



Catalog Number: 116367

## Product name

TSPY3 antibody

## Immunogen

[Human TSPY3 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

TSPY3 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in human testis tissue, LNCaP cells. Observed molecular weight by Western-blot: 39-40 kDa

## Preparation

This antibody was obtained by immunization of TSPY3 recombinant protein (Accession Number: NM\_001077697). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116367(TSPY3 antibody) at dilution of 1:2000



Immunohistochemical of paraffin-embedded human testis using Catalog No:116367(TSPY3 antibody) at dilution of 1:100 (under 10x lens)