

# Anti-Glutaminase Rabbit antibody



Catalog Number: 169077

## Product name

Anti-Glutaminase Rabbit antibody

## Dilutions

WB: 1/1000-1/2000

## Specificity

IHC: 1/50-1/200

Human, Mouse, Rat

ICC/IF: 1/50-1/200

## Antibody description

Glutaminase Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Glutaminase

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal


## Ig Type

Rabbit IgG


## Applications

WB, IHC-P, ICC/IF

## Validations



Western blot analysis of Glutaminase expression in 293T cell lysate.



Immunohistochemical analysis of paraffin-embedded human brain, using Glutaminase Antibody.