

LCORL antibody



Catalog Number: 112177

Product name

LCORL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LCORL Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, PC-3 cells. Positive IHC detected in human lung. Observed molecular weight by Western-blot: 67kd

Preparation

This antibody was obtained by immunization of LCORL recombinant protein (Accession Number: NM_153686). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112177(LCORL antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung using Catalog No:112177(LCORL antibody) at dilution of 1:100 (under 10x lens)