

GINS3 antibody



Catalog Number: 110977

Product name

GINS3 antibody

IF: 1:20-1:200


Validations

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description


GINS3 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in MCF7 cells, HeLa cells. Observed molecular weight by Western-blot: 25-30 kDa



Preparation

This antibody was obtained by immunization of GINS3 recombinant protein (Accession Number: NM_022770). Purification method: Antigen affinity purified.

WB result from Andre Franz in HeLa cell of



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Catalog No:110977(anti-GINS3).


Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110977(GINS3 antibody) at dilution of 1:1000



Ig Type

Rabbit IgG

Applications

ELISA, WB, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:110977(GINS3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)