

# Anti-HCE Rabbit antibody



Catalog Number: 168977

## Product name

Anti-HCE Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

HCE Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human HCE

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP


## Dilutions

WB: 1/1000

ICC/IF: 1/50

IP: 1/20

## Validations



Western blot detection of HCE in K562,C6,3T3 cell lysates using HCE Rabbit mAb(1:1000 diluted).Predicted band size:69kDa.Observed band size:69kDa.