

# Anti-eIF4E Rabbit antibody



Catalog Number: 167277

## Product name

Anti-eIF4E Rabbit antibody

WB: 1/5000

IHC: 1/20

## Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20

## Antibody description

eIF4E Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human eIF4E

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, FC, IP

## Dilutions

## Validations



Western blot detection of eIF4E in K562, Rat Brain, C6, 3T3, HeLa cell lysates using eIF4E Rabbit pAb(1:1000 diluted). Predicted band size: 25kDa. Observed band size: 25kDa.



Immunohistochemistry of eIF4E in paraffin-embedded Human lung cancer tissue using eIF4E Rabbit pAb at dilution 1/20