

BBS9 antibody



Catalog Number: 108377

Product name

BBS9 antibody

WB: 1:200-1:1000

IP: 1:200-1:2000

Specificity

Human, Mouse; other species not tested.

IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

BBS9 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, human heart tissue. Positive IF detected in hTERT-RPE1 cells. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of BBS9 recombinant protein (Accession Number: BC032715). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, IF, WB, IP


Dilutions

Recommended Dilution:

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108377(BBS9 antibody) at dilution of 1:300




Immunohistochemical of paraffin-embedded human liver using Catalog No:108377(BBS9 antibody) at dilution of 1:50 (under 10x lens)


BBS9 antibody




Catalog Number: 108377



Immunohistochemical of paraffin-embedded human liver using Catalog No:108377(BBS9 antibody) at dilution of 1:50 (under 40x lens)



IF result (cytoplasm stain) of anti-BBS9 (Catalog No:108377; 1:50) with hTERT-RPE1 cell by Dr. Moshe Kim.



IP Result of anti-BBS9 (IP:Catalog No:108377, 4ug; Detection:Catalog No:108377 1:500) with mouse testis tissue lysate 4000ug.