

Anti-p60 katanin Rabbit antibody



Catalog Number: 168687

Product name

Anti-p60 katanin Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, ICC/IF

Antibody description

p60 katanin Rabbit monoclonal antibody

Dilutions

1:1000-1:5000

Preparation


Antigen: A synthetic peptide of human p60 katanin

1:20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of p60 katanin in K562, HeLa cell lysates using p60 katanin Rabbit mAb(1:1000 diluted). Predicted band size: 56kDa. Observed band size: 56kDa.