

Anti-Amyloid Precursor Protein Rabbit antibody



Catalog Number: 167387

Product name

Anti-Amyloid Precursor Protein Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Amyloid Precursor Protein Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Amyloid Precursor Protein

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000-1/5000

IHC: 1/20


ICC/IF: 1/20

FC: 1/20


Validations



Western blot detection of Amyloid Precursor Protein in Rat Brain lysates using Amyloid Precursor Protein Rabbit pAb(1:500 diluted). Predicted band size: 87kDa. Observed band size: 100kDa.



Immunocytochemistry of Amyloid Precursor Protein(green) in Hela cells using Amyloid Precursor Protein Rabbit pAb at dilution 1/50, and



DAPI(blue)

Immunohistochemistry of Amyloid Precursor Protein in paraffin-embedded Human breast cancer tissue using Amyloid Precursor Protein Rabbit pAb at dilution 1/100