

Cytohesin 1 antibody



Catalog Number: 109787

Product name

Cytohesin 1 antibody

Immunogen

[Human CYTH1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Cytohesin 1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of Cytohesin 1 recombinant protein (Accession Number: NM_004762). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109787(CYTH1 antibody) at dilution of 1:500