

GRB10 antibody



Catalog Number: 111197

Product name

GRB10 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Antibody description

GRB10 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human kidney tissue. Positive IF detected in HepG2 cells. Positive WB detected in mouse liver tissue, HeLa cells, HepG2 cells, mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 60-70 kDa

Preparation

This antibody was obtained by immunization of GRB10 recombinant protein (Accession Number: XM_017012064). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IHC


Dilutions

Recommended Dilution:


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111197(GRB10 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:111197(GRB10 Antibody) at dilution of 1:50



GRB10 antibody



Catalog Number: 111197

Immunofluorescent analysis of HepG2 cells, using GRB10 antibody Catalog No:111197 at 1:50

dilution and Rhodamine-labeled goat anti-rabbit IgG (red).