

AMOTL1 antibody



Catalog Number: 106997

Product name

Polyclonal

AMOTL1 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

AMOTL1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 107kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_130847). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108000(AMOTL1 antibody) at dilution of 1:500