

ZP2 Polyclonal Antibody



Catalog Number: 160297

Product name

ZP2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ZP2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 651-745 of human ZP2 (NP_003451.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:1000 - 1:4000


IF 1:50 - 1:200

Validations




Western blot - ZP2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ZP2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunofluorescence - ZP2 Polyclonal Antibody

Immunofluorescence analysis of rat oophoroma cells using ZP2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




ZP2 Polyclonal Antibody

Catalog Number: 160297

Immunofluorescence - ZP2 Polyclonal Antibody


Immunofluorescence analysis of rat oophoroma cells using ZP2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



DAPI for nuclear staining.


Immunofluorescence - ZP2 Polyclonal Antibody

Immunofluorescence analysis of mouse oophoroma cells using ZP2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence - ZP2 Polyclonal Antibody

Immunofluorescence analysis of mouse oophoroma cells using ZP2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence - ZP2 Polyclonal Antibody

Immunofluorescence analysis of mouse oophoroma cells using ZP2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence - ZP2 Polyclonal Antibody

Immunofluorescence analysis of mouse oophoroma cells using ZP2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence - ZP2 Polyclonal Antibody

Immunofluorescence analysis of mouse oophoroma cells using ZP2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

ZP2 Polyclonal Antibody



Catalog Number: 160297



Immunofluorescence - ZP2 Polyclonal Antibody

Immunofluorescence analysis of mouse oophoroma cells using ZP2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.