

# FOXG1 antibody



Catalog Number: 110697

## Product name

FOXG1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FOXG1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of FOXG1 recombinant protein (Accession Number: BC035020). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110697(FOXG1 antibody) at dilution of 1:500